
AMOS Manual Breast Pump VARYFLO™

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AMOS 手动吸奶器 VARYFLO™

USER INSTRUCTIONS

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Breast-feeding is the best for babies. You have made the right choice by choosing to breastfeed your new baby. Your milk is perfectly adapted to your baby's needs and contains all the antibodies that protect him/her from certain illnesses and allergies. Breastfeeding is also a quality moment of intense complicity between you and your baby, providing you both with well-being and serenity.

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Please take a few minutes to read these user instructions before using your manual breast pump for the first time.

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1 – When should I use the manual breast pump?

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This manual breast pump has been designed to help you breastfeed in all circumstances:

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Our original Swiss technology combines powerful expression with absolute comfort for the mother. Due to this innovative technique you can adjust the vacuum strength just as you want.

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2 – Presentation of the breast pump parts (*diag. A*)

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3 – Recommendations for use

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As with everything, practice makes perfect. You create suction by activating the pump lever. Start with the lowest vacuum level and raise it by turning the suction strength adjuster 1 until it suits your desires (*diag. B*).
实践能带给您最完美的体验，您可以通过活动此款吸奶器的泵杆来创造出您感到最舒适的吸力。请从最低吸力的档位开始，通过旋转吸力调节器（组件 1）逐渐升高吸力的档位，直到您选择到最适合的档位为止。（见图 B）。

- Choose a calm moment and find a comfortable position.
- 请找到一个平静的环境、舒适的姿势来操作。
- Relax your breasts. This makes it easier to extract the milk.
- 请放松您的乳房，更有助于吸出母乳。
- Pump out the milk from each breast in turn.
- 请轮流从左右两侧乳房吸出母乳。
- Do not worry if, at first, you need to wait a little before the milk flows at a good rate.
- 开始操作时，请不要心急，需要等一小段时间您的母乳才能达到一个很好的流速。

4 – User instructions

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- Wash your hands well and ensure that your breasts are clean.
- 操作前请洗净双手并确保您的乳房也是干净的。
- Also check that the breast pump has been sterilized (*chapter 7*).

- 请检查吸奶器是否已经被消毒（清洁消毒说明请见第 7 章）。
- Place the funnel 5 on the breast in such a way that the nipple is in the middle of the massaging breast shield 6.
- 请将漏斗连接器（组件 5）放置在乳房上，确保乳头在按摩胶垫（组件 6）的中心。
- Ensure that the shield 6 adheres well to your breast so no air is sucked in on the sides.
- 请确保按摩胶垫（组件 6）紧紧贴附在您的乳房上，没有空气可以从边缘吸入。
- Take care to keep the bottle in a vertical position so that the milk flows into it and not into the pump.
- 请小心保持储奶瓶于垂直的位置，从而可以让母乳顺畅得流入储奶瓶而不是吸泵中。
- Begin to pump out your milk adjusting the suction strength adjuster 1 to the lowest level; you need to be gentle with your breasts at first. You can gradually increase the strength according to your physiological needs and your own sensibility (diag. B).
- 在吸奶操作前，请将吸力调节器（组件 1）调节至吸力最低的档位，最初的时候您需要温柔得对待您的乳房。操作过程中，您可以根据您的身体状况以及自身的敏感度逐步得提高吸力。（见图 B）

5 – Storing breast milk

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- Breast milk can be stored up to 3 days (maximum 4°C; do not put it in the door compartment) in the refrigerator using AMOS baby bottles. You can also transfer it into a sterilized container suitable for deep-freezing and place it for a maximum of 6 months (-18°C or colder) in the deep freeze.
- 使用 AMOS 奶瓶或储奶瓶，母乳在冰箱冷藏室最多可保存 3 天（4°C 以下，切不可将母乳放置于冰箱门上或靠近冰箱门的位置）。也可以将母乳转存至消毒过的适用于深度冷冻的储奶容器，在冰箱冷冻室最多可保存 6 个月（-18°C 以下）
- You can store defrosted breast milk in the refrigerator for up to 10 hours. Do not put it back into the deep freezer!
- 解冻以后的母乳在冰箱冷藏室最多可保存 10 个小时。切不可将已解冻过的母乳再放回冰箱冷冻室再次冷冻！
- Breast milk must not be thawed in a microwave oven since you must protect its full nutritional value. Place it in the refrigerator or, if you need it urgently, in warm water (maximum 37°C).
- 请不要使用微波炉对母乳进行解冻或加热，否则将会造成母乳中营养物质的损害或流失。请将母乳放置于冰箱冷藏室慢慢自行解冻，如遇紧急情况急需使用，请使用的温水（最高 37°C）进行解冻及加热。

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CAUTION – Keep all small parts away from children.

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7 – Cleaning and sterilizing

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After each use, the breast pump must be cleaned and all parts in contact with milk must be boiled or sterilized (*diag. A*). This is important to ensure that there is no contamination caused by residual milk.

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8 – Putting the manual breast pump back together

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- 1 Insert the anti-overflow cap 4 correctly into the lid 3 (*diag. C*).
- 1 将防溢帽 (组件 4) 正确得插入吸管盖 (组件 3) 。 (见图 C)
- 2 Push the lid 3 with the mounted anti-overflow cap 4 onto the upper funnel opening 5. Note the orientation of the anti-overflow cap 4 (*diag. D*).
- 2 将带有安装好防溢帽 (组件 4) 的吸管盖 (组件 3) 组装至漏斗连接器 (组件 5) 的上部。过程中请注意防溢帽 (组件 4) 的方向 (见图 D)
- 3 Fit the valve 8 on the valve holder 7 (*diag. E*).
- 3 将阀门 (组件 8) 组装至阀门座 (组件 7) 上。 (见图 E)
- 4 Put the valve holder 7 with the mounted valve 8 into the lower funnel opening 5 (*diag. F*).
- 4 将已组装上阀门 (组件 8) 的阀门座 (组件 7) 组装至漏斗连接器 (组件 5) 的下部。 (见图 F)
- 5 Screw the funnel 5 on the baby bottle 9.
- 5 将漏斗连接器 (组件 5) 拧到储奶瓶 (组件 9) 上。
- 6 Take the pump 2 and insert the aspiration tip of the lid 3 into the pump hole 2 (*diag. G*).
- 6 拿起吸泵 (组件 2) ，将吸管盖 (组件 3) 的吸气嘴插入吸泵 (组件 2) 上对应的吸气孔。 (见图 G)

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BABY BOTTLE 125ml INSTRUCTIONS FOR USE

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




Clean before each use to ensure hygiene. To clean, rinse thoroughly and sterilize all parts before each use by boiling in water for 5 minutes. Completely immerse components in water when boiling. To avoid damage, never clean in a pressure cooker. This product is safe for cleaning in dishwashers – unassembled components in top rack only. Always allow cooling down and remove any remaining water before use. Inspect before each use by pulling the feeding teat in all directions. Throw away at the first signs of damage or weakness. Do not leave a feeding teat in direct sunlight or heat, or leave in disinfectant ("sterilizing solution") for longer than recommended, as this may weaken the teat. Prints / motifs on the surface of the bottle may fade after use due to repeated washing or sterilisation. This is quite normal and completely safe as the inks / colours used for the prints are non-toxic and tested to be safe for food or human contact. Screw teat unit firmly onto bottle when using.

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

For your child's safety and health –WARNING! Always use this product with adult supervision. Keep all components not in use out of the reach of children. Never use feeding teats as a soother. Continuous and prolonged sucking of fluids will cause tooth decay. Always check food temperature before feeding. Take extra care when heating in a microwave oven. Remove all covers and other parts before heating food in a microwave oven. Heating in a microwave oven may produce localised high temperatures. Always stir heated food to ensure even heat distribution and test the temperature before serving. Inspect carefully before each use and discard if damaged. Store in a cool dry place out of reach of children! Never leave children unattended!

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* Recommended typical guidelines. All parts must be completely immersed in water when boiling and not in direct contact with heated surfaces.
Follow instructions provided with your sterilizing and / or warming equipment.

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** Recommended typical guidelines for feeding. Always check temperature of contents before feeding. Please read the enclosed instructions.

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AMOS Europe GmbH